



**SKU:** S9653

## **CURA CORDLESS BED MONITOR WITH DUAL TRANSMITTER CAPABILITY**

**\$599.99**

Simple operation and maintenance.  
No more cords between monitor and pads.  
Can be used with up to six cordless sensor  
pads or other triggering devices.

**PLEASE NOTE: This product needs to be  
combined with a transmitter and  
pressure pad which is sold separately.  
Please see below**

### **Features:**

- Modern ergonomic design
- Monitor up to 7 devices simultaneously  
- a single Premium pad/mat and up to 6  
wireless devices.
- Comes with durable ABS plastic  
enclosure and Caregiver Reset Key.
- Alarm memory and auto re-arm.
- Amber LED visual alarm indicator.
- Caregiver secure reset option.
- Battery operated or optional external  
plug pack.
- Up to 6 months battery life (nominal).
- Easy to access battery compartment.
- Can be connected to a nurse call  
system.

- 3 selectable alarm tones.
- Easy thumbwheel adjustable alarm volume.
- “Pause” feature allows a caregiver to attend to the patient without setting the alarm off.
- Alerts for lost nurse call connection, lost AC power supply & low battery.

### **Safety Features**

- Non resettable “Pad Signal Lost” warning.
- Caregiver secure reset option.
- All safety alerts sound even if alarm volume set to silent, including low battery, AC lost, pad signal lost, nurse call lost & change pad.
- Cable strain relief & breakaway cord technology prevents damage & helps reduce the risk of entanglement when used with Premium Pads & Mats.



Scan QR code to  
view product online

## **ADDITIONAL INFORMATION**

|                              |                                                                           |
|------------------------------|---------------------------------------------------------------------------|
| <b>Brand</b>                 | Cura                                                                      |
| <b>DC operation</b>          | using either a 3 x AA batteries or 6V plug pack                           |
| <b>Dimensions</b>            | 120 x 77 x 40mm                                                           |
| <b>Weight</b>                | 150g without batteries                                                    |
| <b>Range</b>                 | Up to 30m range when used with cordless bed pads and chair pads           |
| <b>Program</b>               | Fully user programmable                                                   |
| <b>Operating temperature</b> | 5 – 40 degrees celcius, relative humidity < 85%                           |
| <b>Alarm tones</b>           | 3 selectable alarm tones                                                  |
| <b>Operating frequency</b>   | 433Mhz fully supervised                                                   |
| <b>Socket</b>                | RJ10 type socket for optional wired bed/chair pad or floor mat connection |
| <b>Indication</b>            | Missing AC/DC power supply detection and indication                       |