



SKU:

HEMECRAFT KINGS CUTLERY HANDLES

~~\$29.99~~ \$29.99

Homecraft Kings Cutlery Handles offer a range of ergonomic grip styles designed for use with Kings modular utensils. They help improve control and comfort for people with limited hand mobility or grip strength.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Type
	S1289	\$24.99		Heavy handle
	S1293	\$19.99		Slim handle
	S1285	\$29.99		Contoured Handle
	S1286	\$19.99		Standard Handle

PRODUCT DESCRIPTION

The Homecraft Kings Cutlery Handles are interchangeable, adaptive grips made for Kings modular cutlery. Available in five distinct designs, each handle style supports a different grip preference or functional requirement – allowing tailored cutlery solutions for independent eating.

These handles are ideal for individuals with arthritis, tremors, hand weakness or limited dexterity. Their varied shapes and weights offer options to improve control, reduce effort, and minimise discomfort during meals.

All handles are designed for compatibility with Kings utensils and are easily cleaned after use.

Available handle types

- **Standard handle** – cylindrical hollow grip for general use
- **Contoured handle** – moulded shape with finger grooves for enhanced grip security
- **Heavy handle** – weighted cylindrical design for those with hand tremors
- **Slim handle** – narrow oval grip; can also be used as a core for built-up sleeves
- **Mushroom handle** – rounded knob-style grip suited for painful or weak joints

Who it's for

- For clients with arthritis, tremor, or weak grip needing custom cutlery support.
- For therapists tailoring dining aids to suit specific grip and motor function needs.

Key features

- **Five ergonomic designs** – accommodate a range of hand conditions and grip preferences
- **Compatible with Kings utensils** – fits securely for seamless use
- **Lightweight or weighted options** – supports either minimal effort or enhanced control
- **Easy to clean** – suitable for daily use in home or care settings
- **Designed for adaptive dining** – promotes comfort, confidence and independence

ADDITIONAL INFORMATION

Standard Handle	110 mm 30 g
Contoured Handle	120 mm 62 g
Small Grip	110 mm 20 g
Large Grip	110 mm 50 g
Oval Grip	100 mm 24 g
Double Ended Handle with Strap	110 mm 36 g
Slim Handle	115 mm 20 g
Mushroom Handle	40 mm 42 g
Custom Made Handle	75 g
Slim Lightweight Handle	110 mm 8 g
Lightweight handle	118 mm 14 g
Lightweight Double Ended Handle with Strap	118 mm 12 g
Heavyweight Handle	110 mm 145 g
Brand	Performance Health