



SKU:

CUTLERY - LIGHTWEIGHT FOAM HANDLED

~~\$99.99~~ \$99.99

These stainless steel foam handled utensils are mounted on plastic stems and sealed into light grey foam handles, to give help for those with weak or painful grip. The foam is comfortable to hold and washable.

Each sold separately



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Product
	S2571	\$29.99		Knife
	S2572	\$29.99		Fork
	S2574	\$29.99		Spoon
	S2576	\$29.99		Teaspoon

Image	SKU	Price	Description	Product
	S2334	\$99.99		Set of 4

PRODUCT DESCRIPTION

ADDITIONAL INFORMATION

Length of Handle	115 mm
Diameter of Handle	28 mm
Knife Weight	55 g
Fork Weight	45 g
Spoon Weight	55 g
Teaspoon Weight	45 g
Brand	ILS ADL