




SKU:

DAY CHAIR BARIATRIC

~~\$1,599.99~~ **\$1,599.99**

The Bariatric Day Chair is a high back orthopaedic chair that features wider seat widths and reinforced frames to suit the heavier client. With adjustable height and seat depth, a contoured backrest and a powder coated steel frame.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Seat Width
	S10234	\$1,199.99		540mm
	S10235	\$1,399.99		590mm
	S10236	\$1,599.99		640mm

PRODUCT DESCRIPTION

The Bariatric Day Chair was created to provide comfort and support to bariatric patients.

ADDITIONAL INFORMATION

Brand	Care Quip
Max SWL	300 kg
Arm Height	250 mm
Back Height	740 mm
Seat Height	430 mm – 580 mm
Seat Depth	490 mm – 550 mm