



**SKU:**

## **DAYS TRI WALKER**

**\$249.99**

The Days Tri Walker provides lightweight, manoeuvrable support for those needing balance assistance. With a three-wheel frame, folding design and included storage bag, it's ideal for indoor and outdoor use.



Scan QR code to  
view product online

### **VARIATIONS**

<b>SKU</b>	<b>Price</b>	<b>Description</b>	<b>Colours</b>
S0233	\$249.99	The Days Tri Walker provides lightweight, manoeuvrable support for those needing balance assistance. With a three-wheel frame, folding design and included storage bag, it's ideal for indoor and outdoor use.	Red
S10137	\$249.99	The Days Tri Walker provides lightweight, manoeuvrable support for those needing balance assistance. With a three-wheel frame, folding design and included storage bag, it's ideal for indoor and outdoor use.	Blue

### **PRODUCT DESCRIPTION**

The Days Tri Walker is a compact, three-wheeled mobility aid designed for individuals who need light walking support and a high level of manoeuvrability. Its narrow frame allows easy steering in tight spaces, making it suitable for use at home or while out and about.

The walker features height-adjustable handles to promote an ergonomic walking posture and loop-style brakes for added safety and control.

Its lightweight aluminium frame folds easily for transport or storage.

A convenient storage bag is included, for carrying personal items securely during errands or daily activities.

### **Who it's for**

- For clients with mild balance concerns who benefit from light, manoeuvrable walking support.
- For therapists prescribing a compact mobility aid with functional features and enhanced turning capability.

### **Key features**

- **Three-wheel frame** – allows nimble movement and tighter turning radius
- **Folding aluminium design** – lightweight and simple to store or carry
- **Height-adjustable handles** – supports safe, ergonomic posture during walking
  - **Loop handbrakes** – provide control and enhance safety
  - **Storage bag included** – helps carry small items while walking

## ADDITIONAL INFORMATION

<b>Colours</b>	Blue, Red
<b>Safe working load (SWL)</b>	125 kg
<b>Gross weight</b>	7.7 kg
<b>Handle height</b>	83 cm – 93 cm
<b>Outer width</b>	68.5 cm
<b>Front castor</b>	7.5"
<b>Rear castor</b>	7.5"
<b>Frame material</b>	Aluminium
<b>Brand</b>	Performance Health

## **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

## **WHAT DIAMETER ARE EACH OF THE WHEELS**

Days 240L Lightweight Aluminium Folding 3 Wheel Tri Walker wheels diameters are 7.5"