



SKU:

DAYS 240L LIGHTWEIGHT ALUMINIUM FOLDING 3 WHEEL TRI WALKER

\$199.99

The Days Aluminium Tri-wheel walker is equipped with a lightweight aluminium frame to ensure users are able to move around easily with the unit. The walker has handles that are easily adjustable in height. The handgrips are perfectly designed for users diagnosed with arthritis. The walker's brakes are easily lockable with a simple push-down technique. The Days lightweight walker is easily foldable using its easy grip handle feature and is suitable for quick storage or transportation.



Scan QR code to
view product online

VARIATIONS

| Image | SKU | Price | Description | Colours |
|--|--------|----------|-------------|---------|
|  A four-wheeled walker with a red frame and a black seat. It has handgrips on the front and a carrying handle on the back. | S0233 | \$199.99 | | Red |
|  A four-wheeled walker with a blue frame and a black seat. It has handgrips on the front and a carrying handle on the back. | S10137 | \$199.99 | | Blue |

PRODUCT DESCRIPTION

The Days Aluminium Tri-wheel Walkers feature a lightweight aluminium frame that is easy to manoeuvre, making them ideal for use indoors or outdoors.

- Ergonomic handgrips for optimum comfort
- Adjustable height handles with 5 settings at 25 mm intervals
- Loop lockable brakes provide added security and convenience. The brake tension can also be easily adjusted to suit requirements.
- Includes durable vinyl bag to transport and conceal your belongings
- Folds easily for storage or transport using the simple easy to grab pull handle
- Optional basket and tray available

ADDITIONAL INFORMATION

| | |
|----------------------------|--|
| Brand | Days Healthcare |
| Maximum user weight | 125 kg |
| Handle height | 830 - 930 mm |
| Overall width | 685 mm |
| Overall depth | 610 mm |
| Unit weight | 5 kg |
| Wheel | 19 cm (7.5") |
| Tyre | EVA |
| Colours | Blue, Red |
| Warranty | 12 Months Warranty Please note: Brake maintenance and adjustment is not covered under warranty. |