



## DAYS HEALTHCARE SERENA ELECTRIC RECLINER LIFT CHAIR - PREMIUM LEATHER - TWIN MOTOR

The Days Healthcare Serena Lift Chair 884 is an electrically operated lift and recline chair which lowers to a fully reclined position. Also available as a Twin Motor chair. Available in Charcoal colour- and made from quality leather material.

### [Days Healthcare Serena Lift Chair Brochure](#)

- Dual-Motor chair features two motors.
- Modern and stylish backrest design.
- It will rise and tilt users to their feet.
- Excellent snooze position for supreme comfort and relaxation.
- Quiet motors provide a smooth and discreet operation.
- Ergonomically contoured seat and back.
- Fully padded chaise for added comfort.
- Quality assured, meeting all relevant fire safety standards.

**EXPRESS DELIVERY TO YOUR DOOR AUSTRALIA WIDE (ONLINE ONLY)**



**SKU:** S2353

**Price:** \$1,699.99 ~~\$2,395.99~~

**Stock:** 15 in stock

### ADDITIONAL INFORMATION

<b>Max SWL</b>	130 kg
<b>Seat Height</b>	500 mm
<b>Seat Depth</b>	520 mm
<b>Seat Width</b>	480 mm
<b>Overall Width</b>	760 mm
<b>Overall Height</b>	1160 mm
<b>brand</b>	Drive Medical
<b>manufacturers-part-number</b>	DAYREC2MBUR





## PRODUCT DESCRIPTION

The Serena Premium Leather 884 Traditional Twin Motor Lift Chair offers great comfort and support with a traditional styling suitable for a wide range of homes. It's dual motor with an early rise ottoman enabling the user to raise their legs without fully reclining the armchair.

### [Days Healthcare Serena Lift Chair Brochure](#)

- Built-in fail-safe battery back-up system.\*
- Comes with a simple-to-operate handset.
- Stitched handset holder to prevent loss / damage to handset.
- Ideal for individuals who have difficulty standing up / sitting down.
- Available in 1 colour option: Charcoal.
- 2-year warranty on wooden frame and electrical components.
- Maximum weight capacity 130kg (20st / 285lb).

## ADDITIONAL INFORMATION

<b>Max SWL</b>	130 kg
<b>Seat Height</b>	500 mm
<b>Seat Depth</b>	520 mm
<b>Seat Width</b>	480 mm
<b>Overall Width</b>	760 mm
<b>Overall Height</b>	1160 mm





**PRICE GUARANTEE** We'll beat it by **5%**  
*If you happen to find a lower price on a stocked item,*

**brand** Drive Medical  
**manufacturers-part-number** DAYREC2MBUR

