


SKU: S2375

## DAYS HEALTHCARE SWIFT WHEELCHAIR

~~\$599.99~~ **\$599.99**






The Days Swift Self-Propelled Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors. The wheelchair is light and is easily folded to allow convenient storage and transport benefits. Available in 3 sizes.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S2785	\$599.99		14"
	S2375	\$549.99		16"
	S1640	\$549.99		18"

## **PRODUCT DESCRIPTION**

### **Key Features**

- Available in 3 Self-Propelled sizes S2785 is 14", S2375 is 16", S1640 is 18"
- Light frame 7.8 kg
- Half-Folding backrest (for storage or transport)
- Quick release wheels (for transport)
- Fixed arms
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Extended brake handle (ease of brake application)
- Upholstery – Black, nylon, padded, fire retardant
- Armrests – Polyurethane foam, fixed, plastic side/skirt guards
- Frame – Aluminium
- Colour – Silver

