



SKU: S2375

DAYS SWIFT SELF-PROPELLED WHEELCHAIR




\$599.99

The Days Swift Self-Propelled Wheelchair is a lightweight, foldable wheelchair ideal for occasional indoor and outdoor use. With three seat width options, it offers ease of transport and storage without compromising basic comfort or functionality.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S2375	\$599.99		16"
	S2785	\$599.99		14"
	S1640	\$599.99		18"

PRODUCT DESCRIPTION

The Days Swift Self-Propelled Wheelchair is a practical solution for short-term mobility support. Its ultra-light aluminium frame makes lifting and handling straightforward, especially during travel or storage.

This model includes a half-folding backrest and quick-release rear wheels for compact storage.

The swing-away, removable leg rests allow for easier transfers, while the solid, puncture-proof tyres provide low-maintenance and reliability on a range of surfaces.

Attendant-controlled drum brakes offer improved control and security, especially when navigating slopes or uneven ground.

With its narrow frame and multiple seat width options, the Days Swift is suitable for tight indoor spaces and general mobility use.

Who it's for

- For clients seeking a lightweight, compact wheelchair for short-term or occasional use.
- For therapists requiring an entry-level mobility aid that's easy to handle, store and maintain.

Key features

- **Ultra-light aluminium frame** – makes lifting and transport simple
- **Half-folding backrest** – supports compact storage
- **Quick-release wheels** – for efficient setup and portability
- **Swing-away, removable leg rests** – assists with safe, efficient transfers
- **Solid, puncture-proof tyres** – reliable and low-maintenance
- **Attendant-controlled drum brakes** – improve braking control on slopes

