



**SKU:** S2070

## **DAYS HEALTHCARE WHIRL WHEELCHAIR**

**\$395.99**

The Days Healthcare Whirl Wheelchair Self-Propelled is designed to be a simple yet versatile wheelchair ideal for part or full-time use. Days Healthcare Whirl's compact frame and dimensions ensure great manoeuvrability indoors and in confined spaces. Its puncture-proof tyres, provide extra peace of mind outdoors.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

The Days Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about.

### **Key Features**

- Quick release wheels - for transport
- Removable armrests
- Swing away and removable leg rests
- Solid tyres
- Upholstery - Nylon
- Colour - Silver

## ADDITIONAL INFORMATION

<b>Brand</b>	Days Healthcare
<b>Seat Width</b>	45 cm
<b>Seat Depth</b>	43 cm
<b>Backrest Height</b>	44 cm
<b>Seat Height (Front)</b>	50 cm
<b>Seat Height (Rear)</b>	47.5 cm
<b>Front Castor Wheel</b>	8"
<b>Rear Wheel</b>	24"
<b>Overall Width</b>	64 cm
<b>Weight</b>	16.8 kg
<b>Max SWL/Weight Capacity</b>	115 kg