


SKU: S2070

DAYS HEALTHCARE WHIRL WHEELCHAIR

\$449.99

The Days Healthcare Whirl Wheelchair Self-Propelled is designed to be a simple yet versatile wheelchair ideal for part or full-time use. Days Healthcare Whirl's compact frame and dimensions ensure great manoeuvrability indoors and in confined spaces. Its puncture-proof tyres, provide extra peace of mind outdoors.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Days Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about.

Key Features

- Quick release wheels - for transport
- Removable armrests
- Swing away and removable leg rests
- Solid tyres
- Upholstery - Nylon
- Colour - Silver

ADDITIONAL INFORMATION

Brand	Days Healthcare
Seat Width	45 cm
Seat Depth	43 cm
Backrest Height	44 cm
Seat Height (Front)	50 cm
Seat Height (Rear)	47.5 cm
Front Castor Wheel	8"
Rear Wheel	24"
Overall Width	64 cm
Weight	16.8 kg
Max SWL/Weight Capacity	115 kg