



SKU: S10383

DECPAC MULTIPURPOSE RAMP

~~\$949.99~~ **\$949.99**



The Decpac Ramp is easily portable and lightweight, allowing patients to use the ramp whenever required. The ramps are made using a combination of carbon materials and unidirectional fiberglass that gives the ramp additional strength – minus the extra weight. The multi-panel design is held together with a patented hinge structure that is also connected by webbing. The ramp folds easily into a concertina style and which makes it easy to use and store away for later use. It is also equipped with a non-slip surface, as well as pads at the top and bottom of the ramp for added safety. The ramp does not require assembly.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S10383	\$849.99		87cm - 5.8kg
	S10384	\$949.99		100cm - 7.5kg

PRODUCT DESCRIPTION

Benefits of the Decpac Ramps

- Easily portable for vehicle transportation and quick storage
- Lightweight and strong yet light enough to be easily transported
- Compact with a concertina system that allows the ramp to fold into a storable size
- The non-slip surface is suitable for all weather conditions and provides added safety
- The ramp is ready for immediate use and does not require assembly
- Wider and longer than a personal ramp and is ideal for travel and commercial usage.

Ramp Sizes

Length	Width	Weight
100cm	83.5cm	NA
87cm	83.5cm	5.6kg

ADDITIONAL INFORMATION

Size	87cm – 5.8kg, 100cm – 7.5kg
Ramp width	83.5cm
Max height capacity	14.5 cm and 16.5 cm
Folded width	20 cm
Ramp weight	5.8 kg and 7.5 kg
Safe working load	300kg