




DECPAC MULTIPURPOSE RAMP

\$570.99 - \$620.99

The Decpac Ramp is easily portable and lightweight, allowing patients to use the ramp whenever required. The ramps are made using a combination of carbon materials and unidirectional fiberglass that gives the ramp additional strength - minus the extra weight. The multi-panel design is held together with a patented hinge structure that is also connected by webbing. The ramp folds easily into a concertina style and which makes it easy to use and store away for later use. It is also equipped with a non-slip surface, as well as pads at the top and bottom of the ramp for added safety. The ramp does not require assembly.

Benefits of the Decpac Ramps :


- Easily portable for vehicle transportation and quick storage
- Lightweight and strong yet light enough to be easily transported
- Compact with a concertina system that allows the ramp to fold into a storable size
- The non-slip surface is suitable for all weather conditions and proves added safety
- The ramp is ready for immediate use and does not require assembly
- Wider and longer than a personal ramp and is ideal for travel and commercial usage.



SKU: S10383

Category: [Wheelchair Ramps](#)

VARIATIONS

Image	SKU	Price	Description	Size
	S10383	\$570.99		87cm - 5.8kg
	S10384	\$620.99		100cm - 7.5kg

GALLERY IMAGES



ADDITIONAL INFORMATION

Size	87cm - 5.8kg, 100cm - 7.5kg
Width	83.5cm
Max Height Capacity	14.5 cm and 16.5 cm
Ramp Weight	5.8 kg and 7.5 kg
Folded Width	20cm
Max SWL	300kg