




SKU:

HERO MEDICAL DELUXE SEAT WALKER WITH ADJUSTABLE HEIGHT

\$199.99

The Hero Medical Deluxe Seat Walker with Adjustable Height has been designed with an easy height adjustment system for the seat and handles, giving users more choice, comfort and safety. It comes with a “light touch” brake system, comfortable padded seat and backrest and can be easily folded for storage and transporting purposes. This Deluxe Seat Walker features a under-seat nylon storage bag, ideal for keeping personal items safe.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
 A four-wheeled walker with a black seat and a blue frame.	S11177	\$199.99		Blue
 A four-wheeled walker with a black seat and a burgundy frame.	S11176	\$199.99		Burgundy

PRODUCT DESCRIPTION

The Hero Medical Deluxe Seat Walker is an excellent walking aid option, available in two colours (Blue and Burgundy). It comes with 6-inch wheels and a 'light touch' braking system which allows users to 'lock' the wheels. The walker's aluminium frame can be adjusted to alter the seat height to maximise comfort, stability and safety. The walker is fitted with ergonomic height adjustable handles to ensure users can maintain an upright posture and to provide ultimate walking support.

Key Features

- Strong lightweight construction
- Height adjustable handles
- Ergonomic handles
- Padded backrest and seat
- User friendly loop brake handles
- "Light touch" braking system
- Simple folding mechanism
- Simple installation
- Nylon storage bag
- Easy manoeuvrability

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight Capacity / Max SWL	140kg (308 lbs)
Overall height (floor to handles)	Min 77cm - Max 101cm
Overall length	Min 62cm Max 76cm
Seat depth	45cm
Seat height	Min 46cm - Max 59cm
Seating width	44cm
Folded width	15cm
Outer width (including back castors)	57cm
Unit weight	6.7kg
Castor size	15 cm/ 6 inch
Warranty	12 months warranty