

SKU: S11788

DIGITENS - TENS & EMS MACHINE

\$199.99

Digital Tens With Soft Case.

Digital Tens and Ems With Soft Bag. Combined EMS and TENS Unit.

Digital Tens has the ability to store your treatment times and schedules which can be examined by your chiropractor. Unit comes with TENS stimulator unit, carrying case, 2X TENS leads, 9V battery, operating manual and packet of 4 Electrodes.



Scan QR code to
view product online



PRODUCT DESCRIPTION

ILS Protens & EMS Dual Unit has fully adjustable parameters and provides relief from acute or chronic pain such as back pain, sports injury, general muscle pain, etc. Battery powered electrical unit which uses electrodes to deliver electrical pulses to the nerve fibres relieving pain.

ADDITIONAL INFORMATION

Brand	ILS physio
Waveform	Asymmetrical Bi-Phasic square pulse
Channel	Dual, isolated between channels
Modes	3 TENS and 3 EMS
Pulse Intensity	Adjustable 0-100mA peak into 500ohm load each channel
Pulse Rate	1 - 150Hz
Pulse Width	30 - 300 micro second
Ramp Time	1 - 30 seconds
Cycle On Time	1 - 30 seconds
Cycle Off Time	1 - 60 seconds
Timer	5 - 95 min. and continuous
Power Source	3 x AAA Batteries (4.5V)