




SKU: S1379

KEY TURNER

\$18.99

This Key Turner is a heavy duty turner designed with a curved built-up handle for comfortable grip to increase leverage while turning a key. The key folds into handle when not in use. The double Key Turner holds two keys and includes a tightening lever that holds the keys in position and simplifies adding or removing keys.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Key Features

- Curved handle for comfortable grip
- Key folds back when not in use
- Holds 2 keys

ADDITIONAL INFORMATION

Brand	ILS ADL
Length	125 mm
Weight	55 grams