




**SKU:**

## **DRESSING STICK**

**\$17.09**

This Dressing Stick assists in the process of dressing and removing garments. Featuring a double wire hook and rubber tip at opposite ends of the stick.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Length
	S1254	\$17.09		18"
	S12745	\$17.09		27"

## **PRODUCT DESCRIPTION**

### Key Features

- Lightweight
- Stylish wooden handle
- Assists in positioning, applying and removing clothing

## ADDITIONAL INFORMATION

<b>Brand</b>	ILS ADL
<b>Length</b>	18", 27"
<b>Weight</b>	87gram, 117gram