



SKU: S12676

BARIATRIC DRIVE MEDICAL ASTON SHOWER COMMODE TRANSIT

\$1,499.99

The Aston Shower commode and hygiene chair is a multipurpose solution for a patient's toileting and showering requirements. The commode chair is perfect for everyday use and features armrests that flip up for side transfers. The Aston commode chair is engineered as a heavy duty model that is easily cleanable and perfect for both commercial and domestic environments.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Aston Shower commode is an everyday multipurpose solution for a patient's toileting and showering needs. The chair features a backrest that adjusts the height to suit the user and can be effortlessly detached for storage and transportation.

The chair's self-propelling wheels provide a quick release to enable storage and transportation and have solid and puncture proof tires.

The Aston Shower commode features,

- Rust-proof and easy reach aluminum clad durable brakes on self-propelling wheels that allow stable positioning.
- The Polyurethane (PU) seat enabled back and armrests are water resistant, easy to clean and comfortable
- Footplates swipe in and out to enable transfers and are easily removable for easy storage and transportation

ADDITIONAL INFORMATION

Width	66cm (26")
Height	97 - 102cm (38" - 40")
Depth	90cm (35")
Seat Width	56cm (22")
Seat Height	51 - 56cm (20" - 22")
Seat Depth	43cm (17")
Weight Capacity	260kg (40st)
Brand	Drive Medical