



**SKU:** S12676

## **BARIATRIC DRIVE MEDICAL ASTON SHOWER COMMODE TRANSIT**

**\$1,499.99**

The Aston Shower commode and hygiene chair is a multipurpose solution for a patient's toileting and showering requirements. The commode chair is perfect for everyday use and features armrests that flip up for side transfers. The Aston commode chair is engineered as a heavy duty model that is easily cleanable and perfect for both commercial and domestic environments.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

The Aston Shower commode is an everyday multipurpose solution for a patient's toileting and showering needs. The chair features a backrest that adjusts the height to suit the user and can be effortlessly detached for storage and transportation.

The chair's self-propelling wheels provide a quick release to enable storage and transportation and have solid and puncture proof tires.

The Aston Shower commode features,

- Rust-proof and easy reach aluminum clad durable brakes on self-propelling wheels that allow stable positioning.
- The Polyurethane (PU) seat enabled back and armrests are water resistant, easy to clean and comfortable
- Footplates swipe in and out to enable transfers and are easily removable for easy storage and transportation

## ADDITIONAL INFORMATION

<b>Width</b>	66cm (26")
<b>Height</b>	97 - 102cm (38" - 40")
<b>Depth</b>	90cm (35")
<b>Seat Width</b>	56cm (22")
<b>Seat Height</b>	51 - 56cm (20" - 22")
<b>Seat Depth</b>	43cm (17")
<b>Weight Capacity</b>	260kg (40st)
<b>Brand</b>	Drive Medical