

SKU:

EASY GLIDE CHAIR - CARER ASSIST

\$949.99

Experience the innovative castor design of the Oscar Easy Glide Carer Assist dining chair. With this braking feature, caregivers can effortlessly assist users in manoeuvring to and from the dining table. Simply apply pressure to the backrest of the dining chair while seated, releasing the castor brakes and enabling smooth movement towards or away from the table.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Colour
	S7867	\$949.99		Coffee
	S7866	\$949.99		Forest Green
	S7865	\$949.99		Pistachio
	S7864	\$949.99		Redgum

PRODUCT DESCRIPTION

The Oscar Easy Glide features a new castor design. This design allows the user to be assisted by a carer in getting to and from a dining table. While seated in the chair, the carer pushes down on the backrest, releasing the castor brakes and allowing the chair to roll toward or away from the table. Sturdy and durable construction, the chair is designed to withstand heavy use and provide reliable support to users.

- Lightweight aluminium frame
- Castors are weight activated to provide stability when getting in and out of the chair
- Woodgrain finish
- Upholstered seat, backrest, and arm pads
- SWL 120 kg

ADDITIONAL INFORMATION

Seat Height	480 mm
Seat Depth	490 mm
Seat Width	450 mm
Height	830 mm
Depth	460 mm
Width	570 mm
Floor to Top of Arm	715 mm
Back Height	370 mm
Max SWL	120 kg