



SKU: B10034

DRIVE MEDICAL EASY RIDER MOBILITY SCOOTER - 2 X 75AH GEL BATTERIES

From: \$7,299.99

The Drive Medical Easy Rider mobility scooter is built for superb comfort, style and performance. It rides and feels much like a motorcycle with features such as telescopic front motorcycle suspension for optimal handling. An easy-to-operate vehicle with many safety features, it will provide you with many years of reliable service.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Drive Medical Easy Rider mobility scooter has the unique design of a classic motorcycle scooter, creating allure and excitement to the experience of driving. It is packed with many safety features and easy-to-use controls.

Key features include:

- Telescopic front motorcycle suspension for optimal handling
- Three large motorcycle style wheels with modern chrome mud guards
- Comfortable rotating sport style seat with sliding and reclining facility
- Adjustable headrest and armrests
- Angle adjustable handle bars
- Throttle lever with automatic safe-braking system
- Additional hand brake for added security
- Illuminated controls feature a twin dial digital dash which allows you to flick between the speedometer, functions, indicator and reset button
- Trip computer with odometer for measuring distance travelled
- All round lighting including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark
- Twin chrome adjustable rear view mirrors which provide extra visibility
- Large accommodating footwell with anti-slip rubber mat
- Key ignition for security
- Rear bumper provides additional protection
- Anti-tip wheels
- Lockable and removable rear pannier - this is an additional accessory which can be purchased separately
- Large seat back storage pocket
- An anti-roll back system that will prevent the scooter from rolling backwards if you are on an incline
- Parking brake for when your scooter is switched off.

ADDITIONAL INFORMATION

Brand	Drive Medical
Range	48 km * **speed and range may vary depending upon user weight, type and incline of terrain, weather, battery charge and condition, operating speed and general driving situation
Top speed	10km per hour
Max SWL	180 kg (396.8 lbs)
Height	115 cm (45 inches)
Length	165 cm (65 inches)
Width	76 cm (30 inches)
Wheel size	45.72 cm (18 inches)
Ground clearance	12 cm (4.7 inches)
Maximum gradient	10 degrees
Battery type	MK Gel 2 x 12V (75AH)
Warranty	2 years warranty for the scooter. 6 months for batteries.