



**SKU:** B10034

## **DRIVE MEDICAL EASY RIDER MOBILITY SCOOTER - 2 X 75AH GEL BATTERIES**

**From: \$7,299.99**

The Drive Medical Easy Rider mobility scooter is built for superb comfort, style and performance. It rides and feels much like a motorcycle with features such as telescopic front motorcycle suspension for optimal handling. An easy-to-operate vehicle with many safety features, it will provide you with many years of reliable service.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

The Drive Medical Easy Rider mobility scooter has the unique design of a classic motorcycle scooter, creating allure and excitement to the experience of driving. It is packed with many safety features and easy-to-use controls.

Key features include:

- Telescopic front motorcycle suspension for optimal handling
- Three large motorcycle style wheels with modern chrome mud guards
- Comfortable rotating sport style seat with sliding and reclining facility
- Adjustable headrest and armrests
- Angle adjustable handle bars
- Throttle lever with automatic safe-braking system
- Additional hand brake for added security
- Illuminated controls feature a twin dial digital dash which allows you to flick between the speedometer, functions, indicator and reset button
- Trip computer with odometer for measuring distance travelled
- All round lighting including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark
- Twin chrome adjustable rear view mirrors which provide extra visibility
- Large accommodating footwell with anti-slip rubber mat
- Key ignition for security
- Rear bumper provides additional protection
- Anti-tip wheels
- Lockable and removable rear pannier - this is an additional accessory which can be purchased separately
- Large seat back storage pocket
- An anti-roll back system that will prevent the scooter from rolling backwards if you are on an incline
- Parking brake for when your scooter is switched off.

## ADDITIONAL INFORMATION

<b>Brand</b>	Drive Medical
<b>Range</b>	48 km * **speed and range may vary depending upon user weight, type and incline of terrain, weather, battery charge and condition, operating speed and general driving situation
<b>Top speed</b>	10km per hour
<b>Max SWL</b>	180 kg (396.8 lbs)
<b>Height</b>	115 cm (45 inches)
<b>Length</b>	165 cm (65 inches)
<b>Width</b>	76 cm (30 inches)
<b>Wheel size</b>	45.72 cm (18 inches)
<b>Ground clearance</b>	12 cm (4.7 inches)
<b>Maximum gradient</b>	10 degrees
<b>Battery type</b>	MK Gel 2 x 12V (75AH)
<b>Warranty</b>	2 years warranty for the scooter. 6 months for batteries.