




**SKU:** B10001

## **DRIVE MEDICAL EASY RIDER MOBILITY SCOOTER - 2 X 100AH GEL BATTERIES**

**From: \$7,999.99**

The Drive Medical Easy Rider mobility scooter is built for superb comfort, style and performance. It rides and feels much like a motorcycle with features such as telescopic front motorcycle suspension for optimal handling. An easy-to-operate vehicle with many safety features, it will provide you with many years of reliable service.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

The Drive Medical Easy Rider mobility scooter has the unique design of a classic motorcycle scooter, creating allure and excitement to the experience of driving. It is packed with many safety features and easy-to-use controls.

### Key features include:

- Telescopic front motorcycle suspension for optimal handling
- Three large motorcycle style wheels with modern chrome mud guards
- Comfortable rotating sport style seat with sliding and reclining facility
- Adjustable headrest and armrests
- Angle adjustable handle bars
- Throttle lever with automatic safe-braking system
- Additional hand brake for added security
- Illuminated controls feature a twin dial digital dash which allows you to flick between the speedometer, functions, indicator and reset button
- Trip computer with odometer for measuring distance travelled
- All round lighting including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark
- Twin chrome adjustable rear view mirrors which provide extra visibility
- Large accommodating footwell with anti-slip rubber mat
- Key ignition for security
- Rear bumper provides additional protection
- Anti-tip wheels
- Lockable and removable rear pannier - this is an additional accessory which can be purchased separately
- Large seat back storage pocket
- An anti-roll back system that will prevent the scooter from rolling backwards if you are on an incline
- Parking brake for when your scooter is switched off.

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	180 kg (396.8 lbs)
<b>Product Battery Size / Type</b>	MK Gel 2 x 12V (100AH)
<b>Product Range</b>	48 km * *speed and range may vary depending upon user weight, type and incline of terrain, weather, battery charge and condition, operating speed and general driving situation
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Black
<b>Product Warranty</b>	Mechanical 12 months Electrical 12 months Batteries 12 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	81 x 147 x 85
<b>Product Size - in carton(s) (gross weight) in kg</b>	149 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	166 x 76 x 115
<b>Product Size - overall assembled (gross weight) in kg</b>	100 kg (excluding Battery)
<b>Product Size - seat base width in cm</b>	52 cm

<b>Product Size - seat base depth in cm</b>	49.5 cm
<b>Product Size - seat back width in cm</b>	51.5 cm
<b>Product Size - seat back height in (inc headrest if applicable) cm</b>	74 cm
<b>Product Size - seat to floor (include range if height adjustable) in cm</b>	67 cm
<b>Product Size - front tyre type/size in inches/solid or pneumatic</b>	19.5" / Pneumatic
<b>Product Size - rear tyre type/size in inches/solid or pneumatic</b>	19.5" / Pneumatic
<b>Product Size - front wheel size in inches</b>	14"
<b>Product Size - rear wheel size in inches</b>	14"
<b>Product Drive - turning circle in cm ( is in radius)</b>	179 cm
<b>Product Drive - max safe ground clearance in cm</b>	12 cm
<b>Product Drive - max safe incline slope in degree</b>	10°
<b>Product Charger Type</b>	8 AH
<b>Product Motor Size</b>	4 Pole / 1300 watt
<b>Product Lighting</b>	LED

**Product Display**

Digital