




SKU:

DRIVE MEDICAL ENVOY 4 PLUS MOBILITY SCOOTER

\$3,299.99

The Drive Medical Envoy 4 Plus combines a stylish design with brilliant driving performance similar to much bigger models. This all-terrain mobility scooter has many excellent features including all-round suspension, intuitive functions and LED headlights. Also featuring long-lasting batteries, it is perfect as a reliable transport option and capable of a 48km travel range.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
 A four-wheeled mobility scooter with a black seat and backrest, a black steering column with a red accent, and a black basket on the front. The body is primarily black with red accents on the front fenders and wheels.	S13070	\$3,299.99		Red
 A four-wheeled mobility scooter with a black seat and backrest, a black steering column with a blue accent, and a black basket on the front. The body is primarily black with blue accents on the front fenders and wheels.	S13071	\$3,299.99		Blue

PRODUCT DESCRIPTION

Drive Medical: Envoy 4 Plus All Terrain Mobility Scooter

Fast, strong and comfortable, the Drive Medical Envoy 4 Plus is ideal for everyday use. It is available in two colours - blue and red.

Key features

- Top speed of 6.5kph
- All-round suspension
- Has a ground clearance of 5cm
- Pneumatic tyres with large 5" rims
- LED front headlight
- A top speed of up to 15kph
- Is powered by 2 x 40AH batteries
- Maximum travel range of 48km per charge
- Two colours to choose from, Red or Blue.

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	160kg (352.74 lbs)
Product Battery Size / Type	2 x 40 AH Lead Acid
Product Range	48 km
Product Speed (Max)	6 km/hr
Product Colour	Blue, red
Product Warranty	Mechanical 12 months Electrical 12 months Batteries 12 months
Product Size - in carton(s) (WxLxH) in cm	62 x 150 x 65
Product Size - in carton(s) (gross weight) in kg	81 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	60 x 121 x 102
Product Size - overall assembled (gross weight) in kg	94kg (207.2 lbs)
Product Size - seat base width in cm	45 cm
Product Size - seat base depth in cm	44 cm

Product Size - seat back width in cm	45 cm
Product Size - seat back height in cm	36 cm
Product Size - seat to floor in cm	52 cm - 57 cm
Product Size - front/ rear tyre type/size	10" / Pneumatic
Product Size - front/rear wheel size in inches	5"
Product Drive - turning circle in cm	135 cm
Product Drive - max safe ground clearance in cm	6.5 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	Yes (Front And Rear CTS Suspension)
Product Charger Type	4 AH
Product Motor Size	2 Pole / 400 Watt
Product Lighting	LED
Product Display	Analog