




SKU:

## **DRIVE MEDICAL ENVOY 4 PLUS MOBILITY SCOOTER**

**\$3,299.99**

The Drive Medical Envoy 4 Plus combines a stylish design with brilliant driving performance similar to much bigger models. This all-terrain mobility scooter has many excellent features including all-round suspension, intuitive functions and LED headlights. Also featuring long-lasting batteries, it is perfect as a reliable transport option and capable of a 48km travel range.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Colour
	S13070	\$3,299.99		Red
	S13071	\$3,299.99		Blue

## **PRODUCT DESCRIPTION**

### **Drive Medical: Envoy 4 Plus All Terrain Mobility Scooter**

Fast, strong and comfortable, the Drive Medical Envoy 4 Plus is ideal for everyday use. It is available in two colours - blue and red.

#### **Key features**

- Top speed of 6.5kph
- All-round suspension
- Has a ground clearance of 5cm
- Pneumatic tyres with large 5" rims
- LED front headlight
- A top speed of up to 15kph
- Is powered by 2 x 40AH batteries
- Maximum travel range of 48km per charge
- Two colours to choose from, Red or Blue.

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	160kg (352.74 lbs)
<b>Product Battery Size / Type</b>	2 x 40 AH Lead Acid
<b>Product Range</b>	48 km
<b>Product Speed (Max)</b>	6 km/hr
<b>Product Colour</b>	Blue, red
<b>Product Warranty</b>	Mechanical 12 months Electrical 12 months Batteries 12 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	62 x 150 x 65
<b>Product Size - in carton(s) (gross weight) in kg</b>	81 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	60 x 121 x 102
<b>Product Size - overall assembled (gross weight) in kg</b>	94kg (207.2 lbs)
<b>Product Size - seat base width in cm</b>	45 cm
<b>Product Size - seat base depth in cm</b>	44 cm

<b>Product Size - seat back width in cm</b>	45 cm
<b>Product Size - seat back height in cm</b>	36 cm
<b>Product Size - seat to floor in cm</b>	52 cm - 57 cm
<b>Product Size - front/ rear tyre type/size</b>	10" / Pneumatic
<b>Product Size - front/rear wheel size in inches</b>	5"
<b>Product Drive - turning circle in cm</b>	135 cm
<b>Product Drive - max safe ground clearance in cm</b>	6.5 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front And Rear CTS Suspension)
<b>Product Charger Type</b>	4 AH
<b>Product Motor Size</b>	2 Pole / 400 Watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog