




SKU:

DRIVE MEDICAL ENVOY 4 PLUS MOBILITY SCOOTER

\$3,299.99

The Drive Medical Envoy 4 Plus Mobility Scooter combines durability with user-friendly functionality. With a 48 km travel range, all-round suspension and intuitive controls, it supports safe and comfortable outdoor mobility.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S13070	\$3,299.99		Red
	S13071	\$3,299.99		Blue

PRODUCT DESCRIPTION

The Drive Medical Envoy 4 Plus Mobility Scooter is a compact yet powerful mobility solution ideal for regular use. Offering the functionality of larger scooters in a streamlined design, it is well-suited to navigating both urban and suburban environments.

Featuring all-round suspension and pneumatic tyres, the Envoy 4 Plus delivers a smoother ride across varied terrain. Its adjustable tiller, armrests and padded seat enable personalised comfort, while delta handlebars enhance steering control.

This scooter supports travel of up to 48 km on a single charge, making it a dependable choice for extended trips. Its LED headlight increases visibility in lower light.

The front basket offers convenient storage for daily essentials.

Available in red or blue, the Envoy 4 Plus offers style alongside its practical benefits.

Who it's for

- For clients seeking a dependable scooter for longer trips with suspension and adjustability.
- For therapists wanting a compact scooter that offers robust range, comfort and safety features.

Key features

- **Up to 48 km travel range** – supports extended daily use without frequent recharging
- **All-round suspension with pneumatic tyres** – ensures smoother rides across uneven surfaces
- **Delta handlebars with tiller adjustment** – improves comfort and manoeuvrability
- **LED front headlight** – enhances visibility for safer travel
- **Padded seat with armrest and height adjustment** – customisable for individual comfort
- **Front basket** – offers practical on-the-go storage
- **160 kg weight capacity** – accommodates a broad range of users
- **Stylish wheels with non-marking tyres** – supports indoor and outdoor use

ADDITIONAL INFORMATION

Brand

Drive Medical