


SKU:

EQUAGEL PROTECTOR CUSHION

~~\$669.99~~ **\$669.99**

The EquaGel Protector Cushion is designed for enhanced postural stability and consistent pressure relief for those at medium risk of injury. Its contoured pelvic zone and temperature-neutral gel construction make it a good choice for those who need reliable positioning and airflow.



Scan QR code to
view product online

VARIATIONS

| Image | SKU | Price | Description | Sizes |
|--|--------|----------|-------------|-----------|
|  | S9967 | \$569.99 | | 16" x 16" |
|  | S9968 | \$569.99 | | 16" x 18" |
|  | S2955 | \$569.99 | | 18" x 16" |
|  | S10283 | \$569.99 | | 18" x 18" |

| Image | SKU | Price | Description | Sizes |
|--|-------|----------|-------------|-----------|
|  | S9970 | \$569.99 | | 20" x 16" |
|  | S9971 | \$569.99 | | 20" x 18" |
|  | S9985 | \$669.99 | | 22" x 16" |
|  | S9986 | \$669.99 | | 22" x18" |

PRODUCT DESCRIPTION

The EquaGel Protector Cushion is a therapeutic seating solution made from 100% EquaGel. Its dual-layer gel structure includes a softer top layer that cushions the body and a firmer base layer that supports posture, helping reduce pressure on sensitive areas.

The cushion provides advanced postural stability, with a 6.4 cm thick EquaGel core with a 1.9 cm contoured pelvic area, which distributes weight evenly during extended sitting. Its design allows for airflow, keeping the seated area cool and dry, and reducing the risk of moisture build-up.

A fluid-resistant, vapour-permeable, flame-resistant cover supports hygiene and safety, while latex-free materials make it suitable for individuals with latex sensitivities.

The EquaGel Protector Cushion also supports a safe working load of 270 kg, accommodating a wide range of individuals.

Who it's for

- For clients at medium risk of pressure injury who need a durable, supportive seating solution.
- For therapists supporting clients who require targeted offloading and enhanced postural stability.

Key features

- **6.4 cm dual-layer gel core** – top layer cushions, while the firm base supports posture
- **Contoured pelvic area** – redistributes pressure and maintains 1.9 cm of gel under key pressure points
- **Fluid-resistant, vapour-permeable, flame-resistant cover** – supports hygiene and user safety
- **Temperature-neutral design** – promotes airflow to keep the seated area cool and dry
- **Latex-free materials** – suitable for those with latex sensitivities
- **Safe working load of 270 kg** – accommodates a wide range of individuals with complex needs

ADDITIONAL INFORMATION

Brand

EquaGel