



SKU:

DAYS FOREARM ERGONOMIC CRUTCHES

\$69.99

The Days Forearm Ergonomic Grip Crutches are lightweight, adjustable crutches designed to enhance user comfort and safety. Featuring ergonomic handgrips and a durable anodised aluminium frame, they are ideal for individuals requiring reliable support.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S1131	\$69.99		Small / Medium
	S1132	\$69.99		Large

PRODUCT DESCRIPTION

Developed in collaboration with physiotherapists, the ergonomic handgrips distribute pressure evenly across the hand, reducing strain during extended use. The closed fixed armband allows the crutch to remain on the forearm when not actively gripped, preventing accidental drops.

Constructed from lightweight anodised aluminium, these crutches offer double adjustability—both the handle-to-armband and handle-to-ground heights can be customised to suit individual user needs.

The non-marking, slip-resistant rubber tips provide added safety on various surfaces.

Who its for

- For clients who require the additional support of a walking aid.
- For therapists prescribing walking aids for greater independence.

Key Features

- **Ergonomic handgrips** – reduce hand strain and enhance comfort during use
- **Double adjustable height** – customise both handle-to-armband and handle-to-ground measurements
- **Closed fixed armband** – keeps crutch in place when not actively gripped
- **Lightweight anodised aluminium frame** – ensures durability and ease of use
- **Non-marking, slip-resistant rubber tips** – provide stability on various surfaces
- **Designed in collaboration with physiotherapists** – ensures user-centric features

ADDITIONAL INFORMATION

Height Range - Medium 890 mm – 1200 mm

Height Range - Large 890 mm – 1370 mm

Max SWL 130 kg

Brand Days

Manufacturer's part number JAN-125AT