




SKU: S8542

INVACARE ETUDE PLUS HI LO HOSPITAL BED

From: \$3,299.99

The Etude Plus care hospital bed, with its attractive aesthetics, flexible design and easy handling is ideal for community care. All the hallmarks and proven features of the popular Etude, such as sturdiness, ease of handling and durability, are very much present in this new community bed.

Bed rails sold separately



Scan QR code to
view product online

PRODUCT DESCRIPTION

Portability and storage are key benefits of the Etude Plus Hospital Bed . Its ease and simplistic design is aided by the fact that it requires no tools to set up or dismantle. With all the hallmarks and proven features of the popular Etude series, the Etude Plus is ideal for the community and home care environment. Head board is an optional accessory.

ADDITIONAL INFORMATION

Max SWL	180 kg
Overall Length	2190 mm
Overall Width	1020 mm
Mattress Length	2000 mm
Mattress Width	900 mm
High - Low Height	330 mm - 730 mm or 400 mm - 800mm
Brand	Invacare