




SKU: S0041

EURO CHAIR - HIGHBACK

\$699.99

The Euro Chair Highback is an orthopaedic chair which comes with a range of great features including an adjustable height and seat depth to suit individual user requirements and contoured backrest for comfort and lumbar support.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S0041	\$699.99		Champagne
	S9937	\$699.99		Slate

PRODUCT DESCRIPTION

The seat angle can be adjusted to provide optimum hip support and comfortable wide padded armrests assist in standing and transfers. Finished in a high quality powder coated steel frame and foam and vinyl upholstery. Available in a range of fabrics and vinyls.

ADDITIONAL INFORMATION

Brand	Royale Medical
Max SWL	200 kg
Seat Width	500 mm
Seat Depth	490 mm – 550 mm
Seat Height	430 mm – 580 mm
Arm Height	250 mm
Back Height	740 mm
Overall Width	700 mm
Weight	15 kg