



**SKU:**

## **EURO CHAIR - HIGHBACK**



**\$699.99**

The Euro Chair Highback is an orthopaedic chair which comes with a range of great features including an adjustable height and seat depth to suit individual user requirements and contoured backrest for comfort and lumbar support.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Colour
	S0041	\$699.99		Champagne
	S9937	\$699.99		Slate

## **PRODUCT DESCRIPTION**

The seat angle can be adjusted to provide optimum hip support and comfortable wide padded armrests assist in standing and transfers. Finished in a high quality powder coated steel frame and foam and vinyl upholstery. Available in 2 colours: champagne and slate.

## **ADDITIONAL INFORMATION**

<b>Brand</b>	Care Quip
<b>Max SWL</b>	200 kg
<b>Seat Width</b>	500 mm
<b>Seat Depth</b>	490 mm - 550 mm
<b>Seat Height</b>	430 mm - 580 mm
<b>Arm Height</b>	250 mm
<b>Back Height</b>	740 mm
<b>Overall Width</b>	700 mm
<b>Colour</b>	Champagne, Slate
<b>Weight</b>	15 kg