


**SKU:** S12747

## **FEMALE URINAL CYGNET**

**\$19.99**

This female urinal can be autoclavable for sanitising and comes with an integrated handle. A 1 litre capacity with 100ml graduation markings helps ensure an easy overnight rest without the need to constantly rise for relief. Anatomically designed to enable use while either standing, sitting, or lying. Can also be chemically cleaned. Made from durable polypropylene.



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

### Key Features

- Lightweight
- Easy to read graduations on the side for measuring
- Clip on lid to prevent spillage

## ADDITIONAL INFORMATION

<b>Brand</b>	ILS ADL
<b>Capacity</b>	1 Litre
<b>Length</b>	31cm
<b>Width</b>	16cm