




SKU:

FLUX WALKER

The Flux walker by Hoggi promotes upright posture and natural walking movements. It provides a stable base from which your child can explore the world. Refreshing new colours have been added to the range – including bright yellow, vibrant pink, green and purple.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The walker provides several options for adjustments – including forward and backward handgrip movement, as well as height and width, which gives the Flux advantages and ensures a unique fit for each individual. The walker comes standard with universal grips, swivel locks and reverse locks. Other accessories include a sling or fold up seat, and hip and bottom pad.

The four sizes will suit from early walkers to young adults, with the largest size weight rated to 100 kg.

- Available in four sizes
- Max weight 100 kg

Colours



Accessories

- Bumper wheels and swivel lock
- Friction brake
- Driving brake and wheel lock
- Anti-tipper
- Dynamic hip control
- Hand disk
- Forearm supports with hand grips
- Sling seat and sling seat rail
- Fold up seat
- Hip/bottom pad
- Basket
- Carrier

