




SKU:

DAYS PATTERNED FOLDING WALKING STICK

\$44.99

The Days Folding Patterned Walking Stick combines functionality with style, offering a lightweight, height-adjustable walking aid adorned with vibrant patterns. Its foldable design ensures easy storage and portability for daily use.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Pattern
 Leopard	S2733	\$44.99	Autumn Leopard	
 Floral	S2734	\$44.99	Floral	
 Blossom	S2736	\$44.99	Tapestry Blossom	
 Paisley	S2735	\$44.99	Tropical Paisley	

PRODUCT DESCRIPTION

The Days folding walking stick features a contoured T-handle for a comfortable grip and is available in various eye-catching designs, including floral, blossom, paisley and leopard patterns. The bold colours not only add a touch of personal style but also make the stick easily identifiable.

Constructed from durable yet lightweight aluminium, the stick folds effortlessly, making it ideal for travel or storage when not in use. Its compact form ensures it fits easily into bags or compartments, providing convenience for users on the go.

With a height adjustment range of 84 cm to 94 cm, the stick accommodates users of varying heights, promoting appropriate posture and reducing strain during use. The anti-slip rubber tip enhances stability on different surfaces, ensuring safe and secure mobility.

Who it's for

- For clients with reduced balance, strength or endurance when walking.
- For therapists prescribing walking aids for greater independence.

Key Features

- **Contoured T-handle** – provides a comfortable and secure grip
- **Foldable design** – allows for easy storage and portability
- **Height adjustable (84 - 94 cm)** – accommodates users of different heights
- **Lightweight aluminium construction** – ensures durability without added weight
- **Anti-slip rubber tip** – enhances stability on various surfaces
- **Available in multiple patterns** – adds a personal touch to the mobility aid

ADDITIONAL INFORMATION

Height adjustment range	840mm - 960mm
Pattern	Floral, Tapestry Blossom, Tropical Paisley, Autumn Leopard
Max SWL	100kg