



SKU: S1122

FOLDING WALKING FRAME - BALL AND ROPE

\$129.99

This folding walking frame is manufactured from lightweight aluminium and is designed to easily fold by pulling the ball at the front of the frame. Guaranteed never to rust. The folding walking frame is fitted with anti-slip rubber tips on the front and back legs for extra grip and stability. The walker is supplied without wheels or skis. Ski Glides and Wheels are sold as pairs and can be [purchased separately](#).



Scan QR code to
view product online

ADDITIONAL INFORMATION

Brand	Days Healthcare
Max SWL	125 kg
Frame width- outside	440 mm - 640 mm
Frame width- inside	350 mm - 470 mm
Height adjustable	870 mm - 920 mm
Style	Walker, Wheel, Ski