




**SKU:** S1122

## **FOLDING WALKING FRAME - BALL AND ROPE**

**\$129.99**

This folding walking frame is manufactured from lightweight aluminium and is designed to easily fold by pulling the ball at the front of the frame. Guaranteed never to rust. The folding walking frame is fitted with anti-slip rubber tips on the front and back legs for extra grip and stability. The walker is supplied without wheels or skis. Ski Glides and Wheels are sold as pairs and can be [purchased separately](#).



Scan QR code to  
view product online

## ADDITIONAL INFORMATION

<b>Brand</b>	Days Healthcare
<b>Max SWL</b>	125 kg
<b>Frame width- outside</b>	440 mm - 640 mm
<b>Frame width- inside</b>	350 mm - 470 mm
<b>Height adjustable</b>	870 mm - 920 mm
<b>Style</b>	Walker, Wheel, Ski