


# FOOD PREPARATION SYSTEM

This durable food preparation and cutting board is ideal for people who have difficulty gripping or who have the use of only one hand.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

The ABS plastic clamp fixes items such as tins, bowls, and the stainless spikes hold vegetables for peeling. Four non-slip rubber feet prevent the board from moving. Suction pads are also provided so the food preparation board can be used on wet surfaces.

## **ADDITIONAL INFORMATION**

<b>Dimensions</b>	315 mm x 295 mm
<b>Weight</b>	1066 g
<b>Brand</b>	ILS ADL