



**SKU:** S1128

## **FOREARM CRUTCHES HEAVY DUTY**

**\$89.99**

The Bariatric Forearm Crutches are extra long in length and can be cut down to suit individual user requirements. Ideal for the bariatric / heavier user who uses crutches for long periods of time.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

A unique contoured plastic handle evenly distributes pressure across the palm of the hand and double thickness reinforced tubing adds extra strength. Aluminium construction reduces weight.

## **ADDITIONAL INFORMATION**

<b>Brand</b>	Sunrise Medical
<b>Handgrip Height</b>	96cm
<b>Cuff Diameter</b>	10cm
<b>Cuff Handle</b>	22.9cm
<b>Max SWL</b>	190 kg