

SKU:

FORTE ICON PLUS MATTRESS (ILS)

~~\$1,299.99~~ **\$1,299.99**

The Forte Icon Plus Mattress delivers reliable static pressure care for individuals at low to medium risk of injury. Exclusive to ILS, this durable mattress balances comfort and clinical performance across aged care, hospital and home settings.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Sizes
	S16398	\$1,199.99	Single	
	S16399	\$1,199.99	Long Single	
	S16401	\$1,299.99	King Single	
	S16400	\$1,299.99	Long King Single	

PRODUCT DESCRIPTION

The Forte Icon Plus Mattress offers consistent pressure relief and comfort in a range of care environments. It features a triple-layered core, including a premium-grade mid-layer foam and an Airoform upper layer designed to increase immersion and redistribute pressure. Profiled stabilising sides promote safer transfers and support central positioning, while the mattress easily conforms to profiling beds. Finished with the Advantiflex medical-grade cover, the Icon Plus meets infection control standards and provides breathable, moisture-resistant protection.

Who it's for

- For clients who require a comfortable support surface for low to medium pressure injury risk.
- For therapists prescribing a reliable mattress with consistent performance for use in aged care, hospital wards or home care.

Key features

- **Airoform immersion upper layer** - alleviates localised pressure points, provides maximum contact with patient for distribution of interface pressure, high temperature stability to deliver consistent immersion, performance and comfort
- **Profiled stabilising sides** - enables ease of primary patient care, facilitates patient transfers, unique hinge points allow mattress to easily conform to varied bed positions, encourages central patient positioning
- **Premium grade mid-layer foam** - open cell high quality foam that provides exceptional comfort and postural support for repositioning and patient care
- **Triple layered core** - antimicrobial treatment applied to foam to maximise longevity and core microclimate
- **Advantiflex cover** - tried and proven medical grade pressure care cover suited for all healthcare settings, excellent standards compliance, infection control, performance and patient care. Pressure injury prevention is delivered through multi-directional stretch qualities and moisture vapour permeability. Standard features include high frequency sealed and welded seams, soft touch PU upper and heavy-duty PVC base

ADDITIONAL INFORMATION

Safe Working Load (SWL)	40 kg - 165kg
Length	198 cm, 203 cm, 198 cm, 203 cm
Width	88 cm, 88 cm, 105 cm, 105 cm
Depth	15 cm
Brand	Forte Healthcare