




**SKU:**

## **INVACARE COMET ALPINE PLUS MOBILITY SCOOTER**

**\$5,999.99**

The Invacare Comet Alpine Plus Mobility Scooter is a powerful and durable option designed for long-range travel over varied terrain. With advanced suspension, safety features and ergonomic seating, it supports comfortable and confident outdoor mobility.



Scan QR code to  
view product online

## VARIATIONS





Image	SKU	Price	Description	Colours
	B10012	\$5,999.99		Blue
	B10015	\$5,999.99		Copper
	B10011	\$5,999.99		Green
	B10014	\$5,999.99		Red

Image	SKU	Price	Description	Colours
	B10013	\$5,999.99		Silver
	B10010	\$5,999.99		White

## PRODUCT DESCRIPTION

The Invacare Comet Alpine Plus Mobility Scooter delivers enhanced performance for those who need reliability across outdoor environments. Built with a robust frame and powerful motor, it handles rough terrain with ease thanks to large 13" pneumatic tyres and a superior suspension system. Safety is prioritised with features such as anti-tip wheels, automatic speed reduction when cornering, fault warning system, reversing alert, brake lights, and a 2-step lever to prevent unintentional freewheeling. The built-in hand brake and high ground clearance add further security.

The scooter's seat offers leather-look upholstery, backrest angle adjustment and a sliding base for personalised comfort and added legroom. The soft-touch control panel includes proportional lever speed control, speed dial, battery level indicator and buttons for lights, horn and indicators. Front and rear LED lighting, splash guards and directional stability ensure better visibility and a smoother experience.

The Comet Alpine Plus provides up to 58 km of range per charge, supporting long journeys with peace of mind.

All Invacare scooters include 12 months of complimentary 24/7 roadside assistance upon registration with the manufacturer.

### Who it's for

- For clients seeking a long-range, stable scooter suited to rough or varied outdoor terrain.
- For therapists needing a robust scooter with advanced suspension and safety systems for their client.

### Key features

- **Up to 58 km range with 160 kg capacity** – ideal for extended outdoor travel
- **Large 13" pneumatic tyres with 4-point suspension** – enhances comfort on uneven surfaces
- **Automatic speed reduction on bends** – improves handling and safety
- **Ergonomic leather-look seat with adjustments** – supports comfort and legroom
- **Anti-tip wheels and 2-step disengage lever** – for added safety
- **Full LED lighting with indicators and brake lights** – for improved visibility
- **Directional stability and fault warning system** – helps maintain safe, straight-line motion
- **User-friendly controls with proportional speed and light functions** – simple and intuitive operation
- **12-month 24/7 roadside assistance** – free with registration through Invacare

## ADDITIONAL INFORMATION

<b>Brand</b>	Invacare
<b>Model</b>	Comet Alpine
<b>Weight Capacity / Max SWL</b>	160 kg
<b>Empty Weight</b>	143 kg
<b>Total Length</b>	145 cm
<b>Driving Unit Width</b>	66.5 cm – 68.5 cm
<b>Total Height</b>	125.5 cm
<b>Floor to Seat Height</b>	44 cm – 51 cm
<b>Backrest Height</b>	63 cm
<b>Backrest Width</b>	51 cm
<b>Seat Depth</b>	46 cm
<b>Front/Rear Tyre Type</b>	Pneumatic
<b>Front/Rear Tyre Size</b>	13' inch
<b>Ground Clearance</b>	10 cm
<b>Turning Diameter</b>	275 cm
<b>Max Safe Slope</b>	10 degrees
<b>Battery Capacity/Size</b>	2x 12V/75Ah
<b>Motor Capacity</b>	650W/1800W
<b>Charger Type</b>	8AH
<b>Maximum Driving Speed</b>	10Km/h

**Maximum Driving  
Distance/Range**

58Km

**Carton Size (W x L x H)**

73 cm x 174 cm x 74 cm