


SKU:

INVACARE COMET ALPINE PLUS MOBILITY SCOOTER

\$5,999.99



The Invacare Comet Alpine Plus Mobility Scooter is a powerful and durable option designed for long-range travel over varied terrain. With advanced suspension, safety features and ergonomic seating, it supports comfortable and confident outdoor mobility.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
B10012	\$5,999.99		Blue
B10015	\$5,999.99		Copper
B10011	\$5,999.99		Green
B10014	\$5,999.99		Red
B10013	\$5,999.99		Silver
B10010	\$5,999.99		White

GALLERY IMAGES



PRODUCT DESCRIPTION

The Invacare Comet Alpine Plus Mobility Scooter delivers enhanced performance for those who need reliability across outdoor environments. Built with a robust frame and powerful motor, it handles rough terrain with ease thanks to large 13" pneumatic tyres and a superior suspension system.

Safety is prioritised with features such as anti-tip wheels, automatic speed reduction when cornering, fault warning system, reversing alert, brake lights, and a 2-step lever to prevent unintentional freewheeling. The built-in hand brake and high ground clearance add further security.

The scooter's seat offers leather-look upholstery, backrest angle adjustment and a sliding base for personalised comfort and added legroom. The soft-touch control panel includes proportional lever speed control, speed dial, battery level indicator and buttons for lights, horn and indicators.

Front and rear LED lighting, splash guards and directional stability ensure better visibility and a smoother experience.

The Comet Alpine Plus provides up to 58 km of range per charge, supporting long journeys with peace of mind.

All Invacare scooters include 12 months of complimentary 24/7 roadside assistance upon registration with the manufacturer.

Who it's for

- For clients seeking a long-range, stable scooter suited to rough or varied outdoor terrain.
- For therapists needing a robust scooter with advanced suspension and safety systems for their client.

Key features

- **Up to 58 km range with 160 kg capacity** – ideal for extended outdoor travel
- **Large 13" pneumatic tyres with 4-point suspension** – enhances comfort on uneven surfaces
 - **Automatic speed reduction on bends** – improves handling and safety
- **Ergonomic leather-look seat with adjustments** – supports comfort and legroom
 - **Anti-tip wheels and 2-step disengage lever** – for added safety
- **Full LED lighting with indicators and brake lights** – for improved visibility
- **Directional stability and fault warning system** – helps maintain safe, straight-line motion
- **User-friendly controls with proportional speed and light functions** – simple and intuitive operation
- **12-month 24/7 roadside assistance** – free with registration through Invacare

ADDITIONAL INFORMATION

Safe Working Load (SWL)	160 kg
Battery Type	Sealed Lead Acid
Battery Size	2x 12V/75 AH
Range	58 km
Max Speed	10 km/h
Size (W x L x H assembled)	66.5 cm – 68.5 cm x 145 cm x 125.5 cm
Weight	143 kg
Suspension	Front & Rear
Display	Digital (Speed, battery, odometer, temperature, time)
Max Safe Incline Slope	10°
Turning Circle	275 cm
Ground Clearance	10 cm
Brand	Invacare

[Invacare Comet Alpine Plus Brochure](#)

[Invacare Comet Alpine Plus User Manual](#)

[Invacare Scooter Specification Chart](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

**HI, COULD YOU LET ME KNOW ABOUT THE COMET ALPHINE SCOOTER'S BATTERIES?
WHAT TYPE ARE THEY?**

The Comet Alphine Mobility Scooter is equipped with 2 x 75AH Gel Batteries