




SKU:

INVACARE COMET ALPINE PLUS MOBILITY SCOOTER

\$5,999.99

The Invacare Comet Alpine is a powerful and reliable mobility scooter designed for outdoor and indoor adventures. It is a sturdy and robust personal mobility device that offers a comfortable and safe ride. With its advanced features and capabilities, the Comet Alpine is an ideal choice for individuals who want to explore the great outdoors with ease and confidence.



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Colours
	B10012	\$5,999.99		Blue
	B10015	\$5,999.99		Copper
	B10011	\$5,999.99		Green
	B10014	\$5,999.99		Red

Image	SKU	Price	Description	Colours
	B10013	\$5,999.99		Silver
	B10010	\$5,999.99		White

PRODUCT DESCRIPTION

Description

The Invacare Comet Alpine mobility scooter features a durable and robust design, with a weight capacity of up to 160 kg. It is equipped with advanced safety features, including anti-tip wheels, a powerful LED lighting system, and an automatic speed reduction system to ensure a safe and secure ride. Comfort is further enhanced with a quality, 'leather look' and feel upholstery seat and backrest. The Comet Alpine is designed for rough terrains, with its pneumatic tires, front and rear suspension system, and high ground clearance. It has a range of up to 58 km on a single charge, making it ideal for long outdoor adventures.

Key features:

- Powerful motors and a superior 4-point independent suspension system to add extra comfort whilst driving
- Directional stability keeps you moving in a straight line.
- Large front and rear pneumatic wheel tyres (13 inch)
- Multiple control features/user-friendly interface - Proportional lever speed control, speed dial control, battery level indicator and soft touch control panel buttons that operate the horn, head lights, warning lights and indicator lights.
- Superior range of safety features - Fault warning system, reversing beeper, low battery warning, hazard lights, indicator lights beeper, speed reduction button, automatic speed reduction on bends, safety flag, a hand brake, brake light, a 2-step disengaging lever to prevent the scooter from freewheeling.
- Built-in anti-splash guards to protect vital part components
- Adjustable seat and backrest - backrest angle adjustment and the seat can slide forward/backward to adjust the depth from the hand controls and create more leg room.
- Front and rear lights and powerful LED lighting system for visibility
- Roadside assistance - If your scooter breaks down you can access the roadside assistance program offered through Invacare to transport you and the scooter to a nominated address. This is a free service for the first 12 months.
- Weight capacity of up to 160 kg

ADDITIONAL INFORMATION

Brand	Invacare
Model	Comet Alpine
Weight Capacity / Max SWL	160 kg
Empty Weight	143 kg
Total Length	145 cm
Driving Unit Width	66.5 cm - 68.5 cm
Total Height	125.5 cm
Floor to Seat Height	44 cm - 51 cm
Backrest Height	63 cm
Backrest Width	51 cm
Seat Depth	46 cm
Front/Rear Tyre Type	Pneumatic
Front/Rear Tyre Size	13' inch
Ground Clearance	10 cm
Turning Diameter	275 cm
Max Safe Slope	10 degrees
Battery Capacity/Size	2x 12V/75Ah
Motor Capacity	650W/1800W
Charger Type	8AH
Maximum Driving Speed	10Km/h

**Maximum Driving
Distance/Range**

58Km

Carton Size (W x L x H)

73 cm x 174 cm x 74 cm