



SKU:

## PRIDE GOGO ULTRA X MOBILITY SCOOTER

~~\$1,999.99~~ **\$1,999.99**

The new and improved Pride GoGo electric travel mobility scooter features an auto-latching lockup mechanism that permits simple frame separation with only one hand. Front frame-mounted seat post offers maximum stability and it's all-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage.

**FREE DELIVERY TO YOUR DOOR.  
AUSTRALIA WIDE (ONLINE ONLY)**



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Wheels
	S6922	\$1,999.99		3 Wheels
		\$1,995.99		4 Wheels

## **PRODUCT DESCRIPTION**

Modular design for easy serviceability and convenient off-board dual voltage charger can charge your battery pack on-board or off-board. At the heart of Pride's philosophy is the understanding that mobility products empower lives. Through industry-leading, consumer-inspired designs, Pride offers an innovative, distinctive and performance driven comprehensive line of mobility products to help customers lead the most active, independent lifestyles.

## ADDITIONAL INFORMATION

<b>Range</b>	Up to 13 km
<b>Speed</b>	Up to 6 km/h
<b>Length</b>	940 mm
<b>Width</b>	495 mm
<b>Heaviest Part</b>	13 kg
<b>Max SWL</b>	118 kg
<b>Brand</b>	Pride
<b>Manufacturer's part number</b>	GOGOULTX3
<b>Tyre Type</b>	Solid
<b>Battery Type</b>	Sealed Lead Acid
<b>Wheels</b>	3 Wheels, 4 Wheels