



SKU:

PRIDE GOGO ULTRA X MOBILITY SCOOTER

~~\$1,999.99~~ **\$1,999.99**

The new and improved Pride GoGo electric travel mobility scooter features an auto-latching lockup mechanism that permits simple frame separation with only one hand. Front frame-mounted seat post offers maximum stability and it's all-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage.

**FREE DELIVERY TO YOUR DOOR.
AUSTRALIA WIDE (ONLINE ONLY)**



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Wheels
	S6922	\$1,999.99		3 Wheels
		\$1,995.99		4 Wheels

PRODUCT DESCRIPTION

Modular design for easy serviceability and convenient off-board dual voltage charger can charge your battery pack on-board or off-board. At the heart of Pride's philosophy is the understanding that mobility products empower lives. Through industry-leading, consumer-inspired designs, Pride offers an innovative, distinctive and performance driven comprehensive line of mobility products to help customers lead the most active, independent lifestyles.

ADDITIONAL INFORMATION

Range	Up to 13 km
Speed	Up to 6 km/h
Length	940 mm
Width	495 mm
Heaviest Part	13 kg
Max SWL	118 kg
Brand	Pride
Manufacturer's part number	GOGOULTX3
Tyre Type	Solid
Battery Type	Sealed Lead Acid
Wheels	3 Wheels, 4 Wheels