




## HAIR WASHING TRAY FOR BED

**\$110.99**

The Hair Washing Tray enables hair to be washed when a patient is lying in bed and is contoured for user comfort. Waste water drains from spout to suitable container, the Hair Washing Tray is made using white moulded plastic.



---

**SKU: S1211**

**Category: [Personal Grooming](#)**

### ADDITIONAL INFORMATION

<b>Length</b>	755 mm
<b>Width</b>	395 mm
<b>Weight</b>	700 g
<b>Brand</b>	ILS ADL