

SKU: S1391

HANDY BAR - CAR TRANSFER BAR

\$59.99

The Handybar – Car Transfer Bar is a compact tool that supports safer vehicle entry and exit. Designed for portability, it provides a secure handhold without modifying the vehicle, and includes emergency features like a seatbelt cutter and window breaker.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Handybar – Car Transfer Bar is a multi-purpose support tool that helps people get in and out of a vehicle safely and independently. It fits securely into the U-shaped striker on a car door frame, providing a stable handhold while you enter or exit the car.

No vehicle modification is required – simply insert, use and remove. Lightweight and portable, it can be easily stored in the glovebox or door compartment.

The Handybar also includes an integrated seatbelt cutter and window breaker, offering additional protection in emergency situations.

Who it's for

- For clients who need extra support when entering or exiting a vehicle.
- For therapists wanting a simple, portable transfer aid for their client.

Key features

- **Secure U-shaped striker fit** – provides a firm handhold with no vehicle modification needed
- **Supports up to 159 kg** – strong enough for a wide range of users
- **Portable and compact** – easy to store and transport between vehicles
- **Emergency seatbelt cutter** – allows for quick release in urgent situations
- **Built-in window breaker** – adds safety and peace of mind while travelling

ADDITIONAL INFORMATION

Overall Length	23 cm
Overall Width	3.8 cm
Product Weight	340 g
Safe Working Load (SWL)	159 kg
Brand	Handy Bar
Handle Length	210 mm
Handle Width	40 mm