



SKU:


HERO MEDICAL COMFY LITE TRANSIT WHEELCHAIR

\$899.99

The Hero Medical Comfy Lite Wheelchair Range, exclusively available through ILS, is a premium, lightweight wheelchair range with ergonomic design and enhanced safety features.

Style meets comfort in this aesthetically designed wheelchair range that offers ease of handling, transport and storage. The Hero Medical Comfy Lite wheelchairs are available in three sizes: 16", 18" and 19".

Hero Medical Comfy Lite Transit Wheelchair



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S309152	\$899.99		16"
	S309153	\$899.99		18"
	S309154	\$899.99		19"

PRODUCT DESCRIPTION

The Hero Medical Comfy Lite Wheelchair [range](#) features a compact and lightweight design for easy storage and transport. These premium wheelchairs offer user comfort, safety and ergonomic features. Available in three sizes: 16", 18" and 19" - all with a safe working load of 135 kg.

Features:

- **Lightweight:** The lightweight aluminium frame makes it highly portable and easy to handle
- **Compact storage and transportation:** With a foldable frame, half-folding backrest, removable leg rests and wheels, it is easy to store and transport, making it a practical choice for smaller spaces.
- **Ergonomic design:** Splayed arm design, grooved padded arm pads and soft padded back, and seat covers maximises user comfort.
- **Stylish:** Matte black finish frame and two-toned back and seat cover gives users a contemporary and refined appearance.
- **Safety first:** Equipped with a side-buckle release seat belt, lockable attendant / knee brakes, and telescopic, adjustable anti-tippers.

Compliance and certification:

- Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact, and Fatigue Testing.
- Compliant to ISO 7176-19 certification - Crash Testing (***S309152** & **S309153** are crash test approved).

ADDITIONAL INFORMATION

Size	16", 18", 19"
Brand	Hero Medical