




**SKU:**

## **HERO MEDICAL FEATHER LITE WHEELCHAIR**




**\$399.99**

The Hero Medical Feather Lite Wheelchair weighs just 8.4 kg and is ideal for short-term mobility and easy travel. With padded seating, swing-away leg rests and safety features, it offers convenience without compromise.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Colours
	S301725	\$399.99		Blue
	S301724	\$399.99		Silver
	S301726	\$399.99		Red

## PRODUCT DESCRIPTION

Designed for portability, the Hero Medical Feather Lite Wheelchair features an ultra-light aluminium frame and half-folding backrest for compact storage.

The swing-away leg rests include foot straps for secure positioning during transit.

This wheelchair includes a seatbelt and anti-tippers for improved safety, while attendant hand brakes provide controlled stopping and parking.

### Who it's for

- For clients requiring a lightweight, transportable wheelchair for short-term or travel use.
- For therapists seeking an affordable transit wheelchair with safety features including anti-tippers.

### Key features

- **8.4 kg aluminium frame** – easy to carry and store
- **Half-folding backrest** – enables compact storage
- **Padded seat and backrest** – for everyday comfort
- **Swing-away leg rests with foot straps** – supports safe positioning
- **Seatbelt and anti-tippers** – enhances stability and safety
- **Attendant hand brakes** – improve control during movement

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Seat width</b>	46 cm (18.1 inches)
<b>Seat depth</b>	40 cm (15.7 inches)
<b>Seat height (From floor)</b>	47.5 cm (18.7 inches)
<b>Backrest (From floor)</b>	85 cm (33.4 inches)
<b>Backrest (From seat)</b>	40.6 cm (15.9 inches)
<b>Armrest height (From floor)</b>	70.5 cm (27.7 inches)
<b>Front castors</b>	6"
<b>Rear wheels</b>	8"
<b>Net weight</b>	8.4 kg (18.5 lbs)
<b>Safe working load</b>	115 kg
<b>Warranty</b>	12 months