




**SKU:**

## **HERO MEDICAL FEATHER LITE WHEELCHAIR 18"**




**\$399.99**

Lightweight and comfortable. The Hero Medical Feather lite Wheelchair 18 inch seat width is a lightweight wheelchair designed to provide comfort and ease of use. The super light (Net Weight 8.4 kg) foldable wheelchair is easy to carry around. Seat belt and attendant hand brakes enhances the safety.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Colours
	S301725	\$399.99		Blue
	S301724	\$399.99		Silver
	S301726	\$399.99		Red

## **PRODUCT DESCRIPTION**

### **Hero Medical Feather Lite Wheelchair 18"**

#### **Stability**

The Hero Medical wheelchair is designed to be safe and stable as long as the centre of gravity is maintained when in use. Avoid sudden movements which may result in shifting your centre of gravity as this will affect the chair's stability.

#### **Inclined Surface**

Be mindful of your strength, capabilities and limitations before attempting to wheel this wheelchair up or down inclined surfaces including ramps. For correct techniques on ramp negotiation, a qualified healthcare professional may be able to provide product demonstration and/or advice.

#### **Hand Brakes**

This wheelchair comes with attendant handbrakes to provide better control and safety.

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Seat width</b>	46 cm (18.1 inches)
<b>Seat depth</b>	40 cm (15.7 inches)
<b>Seat height (From floor)</b>	47.5 cm (18.7 inches)
<b>Backrest (From floor)</b>	85 cm (33.4 inches)
<b>Backrest (From seat)</b>	40.6 cm (15.9 inches)
<b>Armrest height (From floor)</b>	70.5 cm (27.7 inches)
<b>Front castors</b>	6"
<b>Rear wheels</b>	8"
<b>Net weight</b>	8.4 kg (18.5 lbs)
<b>Safe working load</b>	115 kg
<b>Warranty</b>	12 months