



**SKU:**

## **HERO MEDICAL OUTDOOR LITE WALKER - ALUMINIUM**

**\$499.99**

Exclusive to ILS, the Hero Medical Outdoor Aluminium Lite Walker has been designed to maximise your mobility and independence. It has many premium luxury features to ensure you experience effortless mobility outdoors including a lightweight aluminium frame of 6.4 kg\*, large front wheels, sturdy braking and safety systems and easy-to-operate ergonomic functions.

The Hero Medical Outdoor Lite Walker is the perfect blend of lightness and durability. Boasting a unique curved design, the Hero Medical Outdoor Lite Walker is super stylish and is available in two colours (red and black) and three sizes (low, standard and tall).



Scan QR code to  
view product online

## VARIATIONS





Image	SKU	Price	Description	Height	Colour
	S301678	\$499.99		Tall	Black
	S301675	\$499.99		Tall	Red
	S301674	\$499.99		Standard	Red
	S301677	\$499.99		Standard	Black

Image	SKU	Price	Description	Height	Colour
	S301676	\$499.99		Low	Red
	S301679	\$499.99		Low	Black

## PRODUCT DESCRIPTION

The Hero Medical Outdoor Aluminium Lite Walker offers premium luxury features to ensure you experience effortless mobility outdoors. As a heavy-duty quad seat walker, it has a safe working limit of 136 kg.

Key features:

- Lightweight aluminium frame of 6.4 kg\*
- Large 9.5" front castors for easy maneuvering and stability
- Padded easy-fold rounded adjustable backrest for comfortable seating. There are four levels of adjustment which can be easily made by pressing a clip.
- Light touch hand brake system with an integrated parking brake function
- Super easy to fold and store
- Waterproof removable nylon storage bag that straps easily on the seat walker frame
- Spill proof nylon foam padded seat
- Ergonomic and adjustable handle bars with height adjustment clips
- Comes assembled in the box and just needs to be unfolded

\*Note that the gross weight of the product varies:

Low - 7.8 kg

Standard - 7.9 kg

Tall - 8 kg

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Max safe working load (SWL)</b>	136 kg (299 lbs)
<b>Overall height</b>	Low - 82.5 - 95 cm (32.4 - 37.4 inches) Standard - 84.5 - 97 cm (33.2 - 38.1 inches) Tall - 93.5 - 106 cm (36.8 - 41.7 inches)
<b>Seat height</b>	Low - 50.5 cm (19.8 inches) Standard - 52.5 cm (20.6 inches) Tall - 60.5 cm (23.8 inches)
<b>Seat width</b>	Low - 46 cm (18.1 inches) Standard - 46 cm (18.1 inches) Tall - 46 cm (18.1 inches)
<b>Outer width (including back castors)</b>	Low - 60 cm (23.6 inches) Standard - 60 cm (23.6 inches) Tall - 60 cm (23.6 inches)
<b>Inner width (between back frames)</b>	Low - 46.5 cm (18.3 inches) Standard - 46.5 cm (18.3 inches) Tall - 45.5 cm (18.3 inches)
<b>Product frame weight</b>	Low - 6.4 kg (14.1 lbs) Standard - 6.5 kg (14.3 lbs) Tall - 6.6 kg (14.5 lbs)
<b>Gross weight</b>	Low - 7.8 kg (17.1 lbs) Standard - 7.9 kg (17.4 lbs) Tall - 8 kg (17.6 lbs)
<b>Frame material</b>	Alumimium
<b>Colour options</b>	Black and red
<b>Warranty period</b>	12 months for structural frame